

Display Report

Analysis Info

Analysis Name D:\Data\test\0422\qiaosongsu_2_01_1073.d
Method 1073.m
Sample Name qiaosongsu
Comment Agilent ZORBAX SB-C18 5um 4.6*250mm ; Methanol: Water=90:10

Operator
Instrument
bruker
amazon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	140.0 V	n/a	n/a	Trap Drive	50.6
Accumulation Time	1205 μ s	Averages	5 Spectra	Auto MS/MS	off

